



CIP manual

1. Installing CIP Software
2. Installing Image SXM
3. Hardware
4. Microscope
5. Acquire input
6. Prepare input stack
7. Prepare CIP site
8. Running cip1A
9. Running cip1B
10. Running cip2
11. Running cip4
12. Variations

Installing CIP on your Mac

<http://earth.unibas.ch/micro>

1. Download from the website...
...the package from the internet and unzip.
2. Follow the instructions...
... in the README_INSTALL-from-source-file



README_INSTALL-from-source

You will need a fortran compiler, either gfortran or g77.

For Apple OSX, the people from R (<http://r.research.att.com/tools/>) distribute a package of gfortran for OSX. You need developer command line tools installed. If you are lucky you can find a stripped down package as you only need basic tools and dynamic libraries (e.g. libgfortran.dylib, libgcc_s.dylib). If you are unlucky you have to download some 2-4 GByte from Apple to use < 10 MB of tools/libraries. g77 is no longer included in gcc (>3.4).
For Gnu/Linux use gfortran provided from your distribution.
For Microsoft operation systems Cygwin provides a working gfortran or g77. Have a look at <http://gcc.gnu.org/wiki/GFortranBinaries>

To build and install from source:

For OSX:

1) Unpack the archive

2) Edit Makefile and choose your compiler by un/commenting the appropriate line and read the warning

e.g.

```
#COMP = g77
```

```
COMP = gfortran -std=legacy
```

3) To compile all programs type:

```
$ make
```

4) Install all binaries in /usr/local/bin:

```
$ sudo make install
```

IF you have installed the gfortran provided from the link above, no need to worry any further.

Otherwise /usr/local/bin has to exist. To create it type:

```
$ sudo mkdir -p /usr/local/bin
```

/usr/local/bin has to be in your \$PATH variable. If it isn't, modify your .profile file:

```
$ echo "export PATH=/usr/local/bin:$PATH" >> ~/.profile
```

For users of gfortran, using the option -std=legacy is necessary.

Installing Image SXM / Lazy macros and Fiji on your Mac

<http://earth.unibas.ch/micro>

The image is a composite of three screenshots from a Mac OS environment, illustrating the installation of Image SXM and Fiji.


Top Screenshot: A browser window titled "Rock Deformation & Structure Analysis" is open to the website earth.unibas.ch/micro. The page has a navigation menu with "ROCK DEFORMATION" selected. Under "DOWNLOADS", several categories are listed: CIP, FABRIC, GRAIN SIZE, RANDOM, MONTE CARLO, RIG, and LAZY MACROS. A blue arrow points from the "LAZY MACROS" link to a second browser window.

Middle Screenshot: A browser window titled "www.liv.ac.uk/~sdb/ImageSXM/" displays the "Image SXM v 1.94" download page. The page features the text "Image SXM v 1.94 February 2013 Steve Barrett" and a logo of a book. A blue arrow points from the "LAZY MACROS" link in the top screenshot to this window.

Bottom Screenshot: A file browser window titled "LazyMacros" shows 14 items, each a .TXT file. The files are: Lazy ACF Tiling, Lazy CIP LUTs, Lazy D-map, Lazy Digitize, Lazy EBSD, Lazy Erode Dilate, Lazy Grain Boundaries, Lazy Grain Size Mapping, Lazy Lighting, Lazy LUTs, Lazy Pole, Lazy Prepack, Lazy View and Handle, and Lazy Voronoi contacts. A blue arrow points from the "LAZY MACROS" link in the top screenshot to this window.

Bottom Screenshot: A browser window titled "fiji.sc" displays the Fiji software homepage. The page has a title "Fiji Is Just ImageJ" and social media sharing options. The main text describes Fiji as an image processing package based on ImageJ, Java, and Java 3D. It mentions that Fiji is easy to install and has an automatic update function. A "Download Fiji now" button is visible at the bottom left. A blue arrow points from the "LAZY MACROS" link in the top screenshot to this window.

The Hardware



analyzer
fully rotatable

rotating table with
X-Y stage and tilt stage


fully rotatable assembly of
polarizer & λ plate

Use a high quality polarizing microscope.

Use rotating table with mechanical X-Y stage.
Mount tilt stage onto X-Y stage

Use fully rotatable analyzer (= upper polarizer).

Use fully rotatable lower polarizer. Mount full wave plate (λ plate) in 45° orientation with respect to lower polarizer. Wave plate and polarizer have to remain at 45° with respect to each other during rotation.




fully rotatable analyzer (ZEISS)
 λ plate mounted in 45°


block to hold rotatable
analyzer and red filter
(Geowerk, Basel)

red interference filter (660 nm)
inserted below

Insert interference filter with transmittance at 660 nm below lower polarizer. Use 660 nm to be able to observe the images during acquisition etc. If you use a transmittance at 680 or 700 nm, direct observation of images is nearly impossible (700 nm is approximately the spectral limit of human vision).



lower quarter wave
plate (not inserted) lower polarizer
(inserted)




To achieve a state of circular polarization:

- Insert both polarizers (polarizer and analyzer with orthogonal vibration direction)
- Insert quarter wave plate (at 45°) between analyzer and section (use same slot as normally used for full wave plate)
- Insert quarter wave plate (at 45°) between section and lower polarizer (requires additional holder to be mounted above the holder for the polarizer)



Use a tilt stage
(to be acquired from Geoworks Basel)




Use a monochromatic digital camera

Microscope Settings

Hardware specifications:

- cross table with tilt stage
- graticule in eye piece
- red interference filter 660 nm



1. Preparing the microscope

- 1.1. Make sure that the planes of vibration of the polarizers is 0° (horizontal in the image) and 90° (vertical in the image).
- 1.2. Rotate the microscope table to the 0° orientation.
- 1.3. Make sure that the axes of the mechanical X-Y stage are parallel to 0° and 90° .
- 1.4. Make sure that the tilt axes (i.e.the edges) of the tilt stage are parallel to 0° and 90° .
- 1.5. Make sure that the graticule in the eyepiece is parallel to the 0° or 90° direction
- 1.6. Clean & center microscope as best as you can
- 1.7. Place specimen in tilt stage, focus
- 1.8. Carefully install Köhler illumination

2. Use a monochromatic digital camera

- 2.1. make sure that the horizontal and vertical direction of the captured image coincide with 0° and 90° direction.

3. Select sample site and magnification


- 3.1. Match the condenser to the selected magnification of the objective
- 3.2. Do not stop down the condenser diaphragm too much (else the spatial resolution drops)
- 3.3. Select sample site, make sure structural directions match the 0° and 90° direction of the image and the microscope (see section 1 above)
- 3.4. Orient the microscope table parallel to 0° .
- 3.5. Focus image as best as you can (using the focus help of the image grabber)

4. Overview


In order to prepare color micrographs, three interference filters (460, 560, and 660 nm) can be used. The monochrome images are copied into the 3 color channels of Photoshop RGB and then viewed together. Oftenm, the result is very dark, and a gamma of 2.3 is added (⌘L to invoke LUT operations). Using a transmittance of 700 nm instead of 660 nm for the red channel returns slightly more reddish images.

(left:)


Obviously the lower wave plate (plastic foil) is slightly thicker than the upper one (60 μm quartz plate)



R = 660 nm
G = 560 nm
B = 460 nm
using
fixed lambda plate




R = 700 nm
G = 560 nm
B = 460 nm
using
fixed lambda plate



R = 660 nm
G = 560 nm
B = 460 nm
using
rotating lambda plate



Acquiring the input



file.cirpol

1. Circular polarization



- 1.1. Insert lower and upper polarizer.
- 1.2. Insert both quarter wave plates (above and below section)
- 1.3. Insert IF 660 filter
- 1.4. Use linear setting for exposure - do not overexpose



file.nopol

2. No polarization


- 2.1. Remove lower polarizer and both quarter wave plates.
- 2.2. Leave upper polarizer in (to maintain magnification), leave IF 660 filter in.




file.000
...
file.180

3. Rotation images


- 3.1. Insert rotatable lower polarizer / wave plate assembly.
- 3.2. Determine exposure. Check if selected exposure is also OK at other orientations (rotate section). Adapt exposure to brightest orientation (see profiles below).



exposure OK in all orientations



over exposed between 50° and 90° rotation




- 3.3. Carefully return to the exact starting orientation (0°).
- 3.4. Take 19 rotation images with 10° intervals of **CLW** rotation of (upper and lower) polarizer/waveplate assembly

Remember that the exposure selected here remains constant for all rotations and all tilts.


Acquiring the input



eup.stack slice 1




eup.stack slice 2



eup.stack slice 3



eup.stack slice 4



file.eup

4. Tilt images

Four tilts are used.

For each tilt, 3-4 images are acquired with the focus plane being placed moved through the image such that overlapping strips of focus are created.

For each tilt a stack is created.

Using the DepthOf Focus command, each of the tilt stacks is merged to form one completely focused image.

Remember that the selected exposure remains constant for all rotations and all tilts.

import raw image [I]	Insert slice at start [0]
make ROI there [X]	
crop and scale - NN [1]	crop and scale - bilinear [2]
crop and scale - bicubic [3]	
flip stack horizontal [4]	flip stack vertical [5]
rotate entire stack 180° [6]	
revert ROTs in stack	
rotate 180° stack to reference [R]	re-center rotation 0-180 [K]
rotate 90° stack to reference	re-center rotation 0-90
import Text LUT	
apply camera LUT [Y]	invert stack
flare-correct rotated images [F]	REAL flat field [S]
REAL subtract background [Z]	rolling ball [B]
merge tilts [T]	check averages of tilts [M]
smooth stack from to ... [C]	average stack [A]
Grayscale [G]	invert image [J]
scale to pixel [P]	
export to RAW - 19 rot 2 tilts 1 bg [8]	export to RAW - 19 rot 4 tilts 2 bg [9]
make montage [ö]	

Depth Of Focus

Calculation

Near Limit Pixels

Far Limit Pixels

Options


Weighted Average

Reduce Noise First

Acquiring the input

5. Back ground images


5.1. After capturing the last tilt image, the thin section is removed and the background image is acquired.



213mn_10x-071-bg-tilt.tif
mean=217.486

5.2. The section is introduced again in the flat position and focused. Then the section is removed again and a second background image is taken.

NOTE: background images do not have to be that messy ...
In other words: clean the microscope before use - at least, blow out the dust.




213mn_10x-071-bg-flat.tif
mean=216.482

The difference between the mean grey values 1.004 gray values.

The **bg-flat** is used to correct the background of all images except the tilted ones.

The **bg-tilt** is used to correct the tilted images.

Assembling the Input Stack



I. Registration of rotation and tilt images

A separate stack containing

- 19 rotation images: file.000 file.180
 - 4 tilt images: EUP - SUP - WUP - NUP (clockwise)
- is created.

Fiji is used to perfectly match the rotation and tilt images

(In the course of the Stack Registration, the images are inverted)

Credits:

This work is based on the following paper:

P.Thévenaz, U.E. Ruttimann, M. Unser
 A Pyramid Approach to Subpixel Registration Based on Intensity
 IEEE Transactions on Image Processing
 vol. 7, no. 1, pp. 27-41, January 1998.




This paper is available on-line at
<http://bigwww.epfl.ch/publications/thevenaz9801.html>

Other relevant on-line publications are available at
<http://bigwww.epfl.ch/publications/>

Additional help available at
<http://bigwww.epfl.ch/thevenaz/stackreg/>

Ancillary TurboReg_ plugin available at
<http://bigwww.epfl.ch/thevenaz/turboreg/>

You'll be free to use this software for research purposes, but you should not redistribute it without our consent. In addition, we expect you to include a citation or acknowledgment whenever you present or publish results that are based on it.



In Image SXM again, Lazy Prepstack is used to back-invert the rotation and tilt images.

The cirpol and nopol image are introduced at the start (slice 1 and 2), and registered to the first rotation image (slice 3: file.000) by hand.

The nopol image may be flattened using a rolling ball with 100 pixels radius. A duplicate can be used and converted to a mask (white = mask).


The rotation images are now on slice 3-21.

The tilt images are on slice 22-25.

The background images are placed at the end of the stack (slices 26 and 27).

The stack is cropped: here a final size of 1300 · 1000 is used.

The input stack (TIFF) - the input images (RAW)




The input stack

file.stamacrop (stack matched and cropped)


1. cirpol
2. nopol
3. rot.000
4. rot.010
5. rot.020
6. rot.030
7. rot.040
8. rot.050
9. rot.060
10. rot.070
11. rot.080
12. rot.090
13. rot.100
14. rot.110
15. rot.120
16. rot.130
17. rot.140
18. rot.150
19. rot.160
20. rot.170
21. rot.180
22. eup
23. sup
24. wup
25. nup
26. backFlat
27. backTilt

Using the Lazy Prepstack marco, the slices can be exported to single RAW images → CIP




prepare mask and add to input

Preparing a 'CIP site' on the Mac



Prepare CIP site



Prepare sample folder

```

Terminal — bash — 80x25
gpi-051021:~ rheilbronner$
gpi-051021:~ rheilbronner$ cd /Users/rheilbronner/Desktop/2013CIP/CIP\ site\ 4\
tilts\ 13-02-06/CIPdata/prototype
gpi-051021:prototype rheilbronner$ ls
CIP1A      CIP2      input      prototype.ctrl4  prototype.ctrlB
CIP1B      CIP4      prototype.ctrl2  prototype.ctrlA
gpi-051021:prototype rheilbronner$
gpi-051021:prototype rheilbronner$
gpi-051021:prototype rheilbronner$ ls ../CLUT
CIP-P-spectrum.CLUT  CIP-P-standard.POL  crosspol.CLUT.raw
CIP-P-spectrum.POL  CIP-P-standard.pale.CLUT
CIP-P-standard.CLUT  cirpol.CLUT.raw  qtz20-pale.POL
gpi-051021:prototype rheilbronner$ ls ../calib
axio-micro-660-700.CAL  incA-incP-20-660-700.LUT
incA-incP-15-660-700.LUT  incA-incP-20-lin-660-700.LUT
incA-incP-15-lin-660-700.LUT
gpi-051021:prototype rheilbronner$
    
```

Open Terminal

Set path to folder with data
show folder contents

show content of CLUT

show content of calib

Introductory CIP - the '4-tilt intro' run



Prepare control file file.ctrlA for the '4-tilt-intro' run

```


cip-----INPUT-----
1) Title of problem
sample 'prototype' 4 tilts cirpol pre-corrected
2a) x- and y dimension
1300,1000
2b) no.of tilts,tilt-type (2 tilts:"2,012"=soft-EUP-SUP,"2,134"=stiff-WUP-NUP,etc., 4 tilts:"4,0"=soft,"4,1"=stiff)
4,0
2c) Corrections ROT-bg,TILT-bg,FLARE,cirpol-camera? (0=don't l=do)
1,1,1,0
2d) Reference direction for misorientation (N=0,90 E=90,90)
45,52
3) Rotation images (Number = N$ROT):
   at 10 degrees
input/prototype.010
   at 20 degrees
input/prototype.020
   at 30 degrees
input/prototype.030
.....etc.
.....etc.
   at 150 degrees
input/prototype.150
   at 160 degrees
input/prototype.160
   at 170 degrees
input/prototype.170
   at 180 degrees
input/prototype.180
4a) Tilt images (Number = N$TILT):
   First East up
input/prototype.eup
   Second South up
input/prototype.sup
   Third West up
input/prototype.wup
   Fourth North up
input/prototype.nup
4b) Use CirPOL (=N$INCPOL=2) or AMPLitude (=N$INCAMP=3)
2
5) Circular polarisation (Number = 1):
input/prototype.cirpol
6a) Background for ROT etc. images
input/prototype.backR
6b) Background for TILT images
input/prototype.backT
7) Calibration of program and camera
../calib/axio-micro-660-700.CAL
8) Calibration of inclination
../calib/incA-incP-20-lin-660-700.LUT
9) Stereographic Colour Lookup Table
../clut/CIP-P-standard.CLUT
cip-----CIPResults-----
1) Primary result files
   max value
CIP1A/prototype.MAX
   min value
CIP1A/prototype.MIN
   phase of max
CIP1A/prototype.FMAX
   tilt indicator
CIP1A/prototype.Tindex
   error
CIP1A/prototype.ERR
2) Final result files
   azi
CIP1A/prototype.AZI
   inc from polarization
CIP1A/prototype.INCP
   inc from amplitude
CIP1A/prototype.INCA
   edges from polarization
CIP1A/prototype.EDGP
   edges from amplitude
CIP1A/prototype.EDGA
   absolute misorientation from North (0,90)
CIP1A/prototype.MISN
   absolute misorientation from East (90,90)
CIP1A/prototype.MISE
   absolute misorientation from UP (0,0)
CIP1A/prototype.MISH
   absolute misorientation from ref.direction
CIP1A/prototype.MIS-45-52
3) c-axis orientation image (NON-INTERLEAVED)
CIP1A/prototype.CO1A.raw
4) histogram of azi/dip (5° boxes) as from MENTEX -> INVPIMA
CIP1A/prototype.CPF
  
```

description of sample
 x- y- dimension
 no. of tilts = 4, type = soft
 background corrections for rot and tilt, flare correction, no camera correction for cirpol
 azi / inc of octahedral plane


folder/filename of 10° rot
 folder/filename of 20° rot
 folder/filename of 30° rot
 folder/filename of 150° rot
 folder/filename of 160° rot
 folder/filename of 170° rot
 folder/filename of 180° rot
 folder/filename of EastUp
 folder/filename of SouthUp
 folder/filename of WestUp
 folder/filename of NorthUp



folder/filename of cirpol
 folder/filename of background flat
 folder/filename of background tilted
 path to camera calibration file
 path to inclination calibration
 path to color look-up table



folder/filename for derived maximum gray value
 folder/filename for derived minimum gray value
 folder/filename for angle of maximum GV
 folder/filename tilt index 255 if switch, 0 if no switch
 folder/filename for (curve fit) errors



CLUT

folder/filename for AZIMUTH
 folder/filename for INCLINATION from cirpol
 folder/filename for INCLINATION from amplitude
 folder/filename for orientation gradient from INCP
 folder/filename for orientation gradient from INCA
 folder/filename for misorientation w/r North
 folder/filename for misorientation w/r East
 folder/filename for misorientation w/r Heaven
 folder/filename for misorientation w/r to reference direction (octahedral plane)
 folder/filename for c-axis orientation image
 folder/filename for pole density matrix (used for pole figure)

names that have to be provided

Run cip|A for 4 tilts - raw cirpol - INCP - 0-95%

type:
./cip|A if exe is in folder
cip|A if exe is in /usr/local/bin

```
macbook$ cipla
-----
*cip|a: azi, err, inca, incp, max, min, Tindex, Windex
        maximum image size is      1542288
        re-written for macosx and g77      june-06
        inclinations >90 wrap around      june-08
        inclinations <90 wrap (corrected)  feb-10
        writes wrap index into fmax file   mar-11
        cirpol background sub corrected    aug-12
        does 4 tilts                      feb-13
-----
```

name of control file

```
***** calling control
name of control file >
y.ctrla
no.of tilts and tilt type:      4      0
back1,back2,flaresub,camcorr:  1      1      1      1
../../calib/axio-micro-660-700.CAL
../../calib/incA-incP-20-lin-660-700.LUT
-----
```

```
***** calling readfiles
input/prototype.010
input/prototype.020
input/prototype.030
input/prototype.040
input/prototype.050
input/prototype.060
input/prototype.070
input/prototype.080
input/prototype.090
input/prototype.100
input/prototype.110
input/prototype.120
input/prototype.130
input/prototype.140
input/prototype.150
input/prototype.160
input/prototype.170
input/prototype.180
input/prototype.eup
input/prototype.sup
input/prototype.wup
input/prototype.nup
input/prototype.cirpol
input/prototype.backR
input/prototype.backT
-----
```

```
***** calling viewpix
viewpix: x,y coordinates (end = (0,0)) >
0,0
-----
```

no pixel is viewed before calibration

```
***** calling calibrate
-----
***** calling viewpix after calibration *****
viewpix: x,y coordinates (end = (0,0)) >
0,0
-----
```

no pixel is viewed after calibration

input that has to be provided

Run cip1A for 4 tilts - raw cirpol - INCP - 0-95%

inca:

```
+ 16 *
+ 32 ***
+ 48 *****
+ 64 *****
+ 80 *****
+ 96 *****
+ 112 *****
+ 128 *****
+ 144 *****
+ 160 *****
+ 176 *****
+ 192 *****
+ 208
+ 224
+ 240
+ 256
```

0 %	1 %	2 %	5 %	95 %	98 %	99 %	100 %
-----	-----	-----	-----	------	------	------	-------

```
amplitude image incamp
0-1-2-5-95-98-99-100 % = 0 18 28 45 186 189 191 210
input cut-off values (min,max) >
0,186
```

incp:

```
+ 16 *****
+ 32 *****
+ 48 *****
+ 64 *****
+ 80 *****
+ 96 *****
+ 112 *****
+ 128 *****
+ 144 *****
+ 160 *****
+ 176 *****
+ 192 *****
+ 208
+ 224
+ 240
+ 256
```

0 %	1 %	2 %	5 %	95 %	98 %	99 %	100 %
-----	-----	-----	-----	------	------	------	-------

```
inclination image from cirpol
0-1-2-5-95-98-99-100 % = 0 7 14 28 174 178 181 230
input cut-off values (min,max) >
0,174
```




type GV of 0% and 95% cut-off

```
-----
***** calling writeprimary
-----
***** calling writefiles
-----
***** calling polefigure
-> maximum of polefigure is at (25, 5): 6.26243
-----
***** calling ava
polamp=2 => cirpol used for coi
name of clut:
.././clut/CIP-P-standard.CLUT
```









maximum within half halfspace

input that has to be provided

CIPIA results for 4 tilts - raw cirpol - INCP - 0-95%



CIPIA results for 4 tilts - raw cirpol - INCP - 0-95%



Prepare control file file.ctrlB for the '4-tilt-intro' run

```

cip-----INPUT-----
1) Title of problem
sample 'prototype' 4 tilts cirpol pre-corrected
2a) x- and y dimension
1300,1000
2b) no.of tilts,tilt-type (2 tilts:"2,012"=soft-EUP-SUP,"2,134"=stiff-WUP-NUP,etc., 4 tilts:"4,0"=soft,"4,1"=stiff)
4,0
2c) Corrections ROT-bg,TILT-bg,FLARE,cirpol-camera? (0=don't 1=do)
1,1,1,0
2d) Reference direction for misorientation (N=0,90 E=90,90)
45,52
3) azi image:
CIP1A/prototype.AZI
4) inc0 images (0-90):
CIP1A/prototype.INCP
5) T-index file (switch where GV > 0):
CIP1A/prototype.Tindex
7) Calibration of program and camera
../../calib/axio-micro-660-700.CAL
8) Calibration of inclination
../../calib/incA-incP-20-lin-660-700.LUT
9) Stereographic Colour Lookup Table
../../clut/CIP-P-standard.CLUT
cip-----CIP1results-----
1) Final result files
CIP1B/prototype.INC
2) c-axis orientation image (NON-INTERLEAVED)
CIP1B/prototype.COI1B.raw
3) histogram of azi/dip (5° boxes) as from MENTEX -> INVPIMA
CIP1B/prototype.CPF

```

description of sample
x- y- dimension
no. of tilts = 4, type = soft
background corrections (see prototype.ctrlA)
azi / inc of octahedral plane
folder/filename of azimuth
folder/filename of preliminary inclination 0°-90°
folder/filename of Tilt Index file
path to camera calibration file
path to inclination calibration
path to color look-up table
folder/filename for derived full inclination 0°-180°
folder/filename for c-axis orientation image
folder/filename for pole density matrix (used for pole figure)



names that have to be provided

Run cip1B for 4 tilts - raw cirpol - INCP - 0-95%

type:
./cip1B if exe is in folder
cip1B if exe is in /usr/local/bin

```
macbook$ cip1b
-----
*cip1b: azi/inc from azi/inc/tindex
      no polefigure !
      maximum image size is      1542288
      re-written for macosx and g77      june-06
      last update                      march-07
-----

***** calling control
name of control file >
y.ctrlb
.././calib/axio-micro-660-700.CAL
.././calib/incA-incP-20-lin-660-700.LUT
-----

***** calling readfiles
CIP1A/prototype.AZI
CIP1A/prototype.INCP
CIP1A/prototype.Tindex
-----

***** calling convert2
***** calling writefiles
-----




***** calling ava
.././clut/CIP-P-standard.CLUT
```

name of control file





reading input files

input that has to be provided


CIPIB results for original cirpol - 0-95% - INCP



soft 4 tilts



stiff 4 tilts



Prepare control file file.ctrl2 for the '4-tilt-intro' run

```

cip-----INPUT-----
1) Title of problem
sample 'prototype' 4 tilts cirpol pre-corrected
2a) x- and y dimension
1300,1000
3) Reference direction for misorientation (N=0,90 E=90,90)
45,90
4) Masking for pole figure and orientation image ? (1=yes, 0=no)
1
5) Want misorientation, edge and orientation image ? (1=yes, 0=no)
1,1,1
6) Pole figure correction 1=sin(inc) 2=sin(inc-1/2) 3=-5°-95° 4=sqrt(), 5=sin(delta)
1
7) Stereographic Colour Lookup Table
../clut/CIP-P-spectrum.CLUT
cip-----INPUT FILES-----
1) azimuth file
CIP1A/prototype.azi
2) inclination file
CIP1B/prototype.inc
3) masking file (0 where o.k. >0 where masked) (different from CIP1 !!!)
input/prototype.mask
cip-----OUTPUT IMAGES-----
1) edges using 2 neighbours (different from CIP1 !!!)
CIP2/prototype.EDG2s
2) edges using 4 neighbours (different from CIP1 !!!)
CIP2/prototype.EDG4a
3) absolute misorientation from North (0,90)
CIP2/prototype.MISN2
4) absolute misorientation from East (90,90)
CIP2/prototype.MISE2
5) absolute misorientation from UP (0,0)
CIP2/prototype.MISH2
6) absolute misorientation from ref.direction
CIP2/prototype.MIS-45-90
7) c-axis orientation image (NON-INTERLEAVED)
CIP2/prototype.COI2.raw
8) histogram of azi/dip (5° boxes) as from MENTEX -> INVPIMA
CIP2/prototype.CPF2
9) pole figure image (36*36) -> Lazy Pole
CIP2/prototype.PFIG2

```

description of sample

x- y- dimension

reference orientation

masking yes

want all types of mis-/orientation images

standard pole density correction

path to color look-up table

input: azimuth file

input: inclination file

input: mask

folder/filename for OGI 2 neighbours

folder/filename for OGI 4 neighbours

folder/filename for MOI w/r to North

folder/filename for MOI w/r to East

folder/filename for MOI w/r to Heaven

folder/filename for MOI w/r to reference direction

folder/filename for COI

folder/filename for pole density matrix (180 · 180)

folder/filename for stereo pole figure matrix (36 · 36)



names that have to be provided

Run cip2 for 4 tilts - raw cirpol - INCP - 0-95%

type:
./cip2 if exe is in folder
cip2 if exe is in /usr/local/bin

macbook\$ cip2

```
-----  
cip2: coi/mis/edg from azi/inc/mask  
      maximum image size is      6000000  
      re-written for macosx and g77      apr-04  
>>> reads square (CLUT) and stereo (POL)  
edg2s sum of difference with 2 neighbours (0-180)  
edg4a 2*av of difference with 4 neighbours (0-180)  
      last update (polfig correction)      march-07  
      last update (edges)      august-07  
+MASK considers mask in misor and edge      june-08  
-----
```

```
*calling control  
name of controle file :
```

name of control file

y.ctrl2

```
-----  
*calling readfiles  
*calling readfiles: azi  
CIP1A/prototype.azi  
*calling readfiles: inc  
CIP1B/prototype.inc  
*calling readfiles: mask  
input/prototype.mask  
*   xdim*ydim=itot      1300      1000      1300000  
-----
```

```
*calling segment  
-----
```

```
*calling writefiles  
CIP2/prototype.EDG2s  
CIP2/prototype.EDG4a  
CIP2/prototype.MISN2  
CIP2/prototype.MISE2  
CIP2/prototype.MISH2  
CIP2/prototype.MIS-45-90  
-----
```

```
*calling polefigure using mask  
*   imask =      1 (0=no,1=yes)  
*   polcorr =      1 (option 1-5)  
--> maximum of polefigure is at (25, 5):
```


5.53562

maximum of pole figure
(need to note this number)



```
-----  
*calling ava  
../../clut/CIP-P-spectrum.CLUT
```

input that has to be provided



CIP2 results for original cirpol - 0-95% - INCP




soft 4 tilts




max = 5.5 / 8



stiff 4 tilts



max = 5.8 / 8



CIP2 results for original cirpol - 0-95% - INCP

SOFT - unmasked



edg2s




edg4a



mis-45-90



misE2




misH2




misN2

CIP2 results for original cirpol - 0-95% - INCP


SOFT - masked




edg2s




edg4a




mis-45-90



misE2



misH2



misN2

Prepare control file file.ctrl4 for the '4-tilt-intro' run

```

cip-----INPUT-----
1) Title of problem
sample 'prototype' 4 tilts cirpol pre-corrected
2a) x- and y dimension
1300,1000
3) 4 reference directions for misorientation (N=0,90 E=90,90)
45,52,135,52,45,128,135,128
4) Masking for pole figure and orientation image ? (1=yes, 0=no)
0
5) Want misorientation, edge and orientation image ? (1=yes, 0=no)
1,1,1
6) Pole figure correction 1=sin(inc) 2=sin(inc-1/2) 3=-5°-95° 4=sqrt(), 5=sin(delta)
1
7) Stereographic Colour Lookup Table
../../clut/qtz20-pale.POL
cip-----INPUT FILES-----
1) azimuth file
CIP1A/prototype.AZI
2) inclination file
CIP1B/prototype.INC
3) masking file (0 where o.k. >0 where masked) (different from CIP1 !!!)
input/prototype.mask
cip-----OUTPUT IMAGES-----
1) edges using 8 neighbours average (different from CIP1 !!!)
CIP4/prototype.EDG8a
2) edges using 8 neighbours max (different from CIP1 !!!)
CIP4/prototype.EDG8m
3) absolute misorientation from North (0,90)
CIP4/prototype.MISr1_045_052
4) absolute misorientation from East (90,90)
CIP4/prototype.MISr2_135_052
5) absolute misorientation from UP (0,0)
CIP4/prototype.MISr3_045_128
6) absolute misorientation from ref.direction
CIP4/prototype.MISr4_135_128
7) c-axis orientation image (NON-INTERLEAVED)
CIP4/prototype.COI4.raw
8) histogram of azi/dip (5° boxes) as from MENTEX -> INVPIMA
CIP4/prototype.CPF4
9) pole figure image (36*36) -> Lazy Pole
CIP4/prototype.PFIG4

```

description of sample

x- y- dimension

4 reference orientations

masking no

want all types of mis-/orientation images

standard pole density correction

path to color look-up table (POL form)

input: azimuth file

input: inclination file

input: mask

folder/filename for OGI 8 neighbours

folder/filename for OGI 8 neighbours

folder/filename for MOI w/r to reference direction 1

folder/filename for MOI w/r to reference direction 2

folder/filename for MOI w/r to reference direction 3

folder/filename for MOI w/r to reference direction 4

folder/filename for COI

folder/filename for pole density matrix (180 · 180)

folder/filename for stereo pole figure matrix (36 · 36)



names that have to be provided

Run cip4 for 4 tilts - raw cirpol - INCP - 0-95%

type:
./cip4 if exe is in folder
cip4 if exe is in /usr/local/bin

```
macbook$ cip4
-----
*cip4: coi/mis/edg from azi/inc/mask (=cip2)
      4 reference directions (*.ctrl4)
      maximum image size is 6000000
      maximum image width is 3000
      re-written for macosx and g77 mar-04
  >>> reads square (CLUT) and stereo (POL)
edg8a 2*av of difference with 8 neighbours (0-180)
edg8m max of difference with 8 neighbours
      last update (polefig corr) march-07
      last update (edges) august-07
+MASK considers mask in misor and edge june-08
-----
```

```
*calling control
name of controle file >
y.ctrl4
```

name of control file

```
-----
*calling readfiles
*calling readfiles: azi
CIP1A/prototype.AZI
*calling readfiles: inc
CIP1B/prototype.INC
*calling readfiles: mask
input/prototype.mask
* xdim*ydim=itot 1300 1000 1300000
-----
```

```
*calling segment
-----
```

```
*calling writefiles
CIP4/prototype.EDG8a
CIP4/prototype.EDG8m
CIP4/prototype.MISr1_045_052
CIP4/prototype.MISr2_135_052
CIP4/prototype.MISr3_045_128
CIP4/prototype.MISr4_135_128
-----
```

```
*calling polefigure using mask
* imask = 1 (0=no,1=yes)
* polcorr = 1 (option 1-5)
--> maximum of polefigure is at (25, 5):
```


5.53562

maximum of pole figure
(need to note this number)


```
*calling ava
../../clut/qtz20-pale.CLUT
-----
```

input that has to be provided


CIP4 results for original cirpol - 0-95% - INCP




soft no mask 4 tilts




max = 4.9 / 8




qtz20-gold.POL




soft with mask 4 tilts



max = 5.5 / 8




qtz20-pale.CLUT




CIP4 results for original cirpol - 0-95% - INCP

SOFT - unmasked



edg8a



edg8m



misr1_045_052



misr2_135_052




misr3_045_128



misr4_135_128

Variations on the theme

Pre-correcting the cirpol



Attempting the impossible

As soon as the section is removed from the microscope set up for circular polarization, extinction is obtained: the image is dark.

A 'standard background image for circular polarization is obtained by taking an image of a lambda plate or similar under circular polarization. This would then only show the distribution of light intensity across the field of view.

Images of thin sections under cross polarizers may show intensity variations (above and beyond the local crystal orientation) that is due to
a section thickness and
b uneven illumination



A central dome may therefore be due to central light dome and/or a central 'dome of thin section thickness'. These two cannot be separated.

The best - or most pragmatic correction is therefore one that corrects for both of these at once.


The following procedure is recommended:

1. Apply camera Palette to cirpol image (use Lazy Lighting).
2. Use Rolling Ball with radius = 100 and the option 'Faster' unchecked.

This pre-corrected cirpol will need NO camera correction in cip I.A.






The reasoning behind this approach is that the texture should be the same everywhere in the thin section. In other words, on average, the texture should be the same everywhere. Brighter areas in a thin section are not attributed to a 'nest' of flatter lying axes but to a greater section thickness in that area and / or to brighter illumination at that spot (a true light dome).




DO NOT...

... use the Rolling Ball on the original micrograph before applying the camera LUT





CIP results for pre-corrected cirpol - 0-95% - INCP




soft 4 tilts



max = 6.3 / 8



stiff 4 tilts



max = 7.4 / 8

